

Product Code . DS-P-11147

Galvanometer



Description

Catalogue No. : "IDAC4267/3.7"

DESCRIPTION :-

- Centre Zero Universal Meter Galvanometer
- The centre zero universal meter (Universal Galvanometer) is electrically identical to the universal meter.
- 1000-ohm movement, moving coil, diode protection against overload up to 12V D.C.
- Range: 50uA-0-50uA.
- The unit is ideal for use as an indicator in slow A.C

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/