



Product Code . DS-P-11152

Battery of Eliminator

Description

Catalogue No. : “IDAC6287/4.2”

DESCRIPTION :-

- Power Supply: 1.5-12V/3A
- Battery Eliminator offered by us is ideally designed to be safe reliable source of low voltage D.C., for use in class rooms.
- It allows teachers to conduct experiments requiring a variety of direct voltages while avoiding the decreased characteristics of dry cell.
- It comprises a heavy-duty transformer wound with copper wire with high lamination to avoid losses, selector switch bridge rectifiers with jewel light and safety fuse.
- Precise results
- Long life
- Affordable

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/