

Product Code . DS-P-11154

One Way Key



Description

Catalogue No. : “IDAC1302/4.4”

DESCRIPTION :-

- Consists of two brass blocks with socket terminals mounted on an insulated base.
- Plastic molded
- Contact is made by a tapered plug with an insulating top.
- Base size: 100x62mm

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/