


Product Code . DS-P-11165

Single Turn



Description

Catalogue No. : "IDAC 0102/5.1"

DESCRIPTION :-

- Magnetic Field Demonstration OHP Single-Turn.
- Magnetic field demonstration coil to demonstrate field distribution of current carrying conductor.
- The pattern can be demonstrated on a large screen using an overhead projector.
- A wide range of demonstration coils having different number of turns is available.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/