



Product Code . DS-P-11169

## Compass

### Description

---

Catalogue No. : "IDAC3613/5.5"

#### DESCRIPTION :-

- Magnetic Compass in Aluminum Case
- Size: 16mm
- Top Transparent.
- With needle pivoted between top and bottom glass plate, without dial, diameter 20mm
- Diameter : 18 mm
- Both side Glass

---

Didac Scientific Pvt. Ltd,  
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom  
Direct Contact Details ✉ [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)  
💻 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)