



Product Code . DS-P-11173

Demonstration Electric Motor with Magnet

Description

Catalogue No. : "IDAC0321/09"

DESCRIPTION :-

- Suitable for demonstration of a DC motor, having a two-pole armature wound with enameled copper wire and a detachable permanent bar magnet (NS indicated) mounted on a base.
- External connections via 4mm sockets.
- Operating voltage: 4-6 V DC small and compact
- Height: 125mm

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/