

**Product Code . DS-P-11174**

# Electrostatic Kit

## Description



**Catalogue No. : “IDAC0023/010”**

### DESCRIPTION :-

- The kit is used for demonstrating a number of electrostatics experiments.
- It consists of various elements and a gold leaf electroscope with a transparent front window and ground glass rear window.

SCOPE OF SUPPLY:

- Proof plane on insulating handle Metalized polystyrene sphere
- Electrophorus plate insulating handle Polythene strip
- Polythene tiles
- Nylon reel
- Aluminum cans
- Cellulose acetate strip Wire stirrup
- Polishing cloth Plated wire hook
- Gold leaf electroscope

Didac Scientific Pvt. Ltd,  
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom  
Direct Contact Details  [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)  
 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)