



Product Code . DS-P-11184

Vernier Caliper

Description

Catalogue No. : "IDAC1447/13.1"

DESCRIPTION :-

- A Vernier Caliper is a visual aid to take an accurate measurement reading between two graduation markings on a linear scale by using mechanical interpolation; thereby increasing resolution and reducing measurement uncertainty by using Vernier acuity to reduce human estimation error.
- Provided with thumb movement for easy motion of jaws.
- Material: steel sheet with chrome plating.
- For internal, external & depth measurements up to 150mm x 0.1mm

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/