

Product Code . DS-P-11188

Spherometer



Description

Catalogue No. : "IDAC4162/014.1"

DESCRIPTION :-

- Range: 1/100MM
- It consists of 3 legs.
- Scales and head disc are made of Brass.
- Vertical scale thickness 1mm.
- Diameter of micrometer head disc 38 mm.
- Head disc is made of brass.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/