

Product Code . DS-P-11203

## Beam Balance



### Description

---

Catalogue No. : "IDAC3234/16.1"

#### DESCRIPTION :-

- It consists of a Brass beam engraved with a rider scale.
- Provided with Agate wedge, detachable pans and housed in a wooden case having glass windows.
- Capacity: 200gm.
- Sensitivity: 5mg

---

Didac Scientific Pvt. Ltd,  
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom  
Direct Contact Details ✉ [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)  
💻 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)