



**Product Code . DS-P-11216**

**Force Table SN720**

## Description

---

Catalogue No. : "IDAC7702/019"

DESCRIPTION :-

- To demonstrate the vector addition of forces.
- Dual scale allows students to take reading clockwise or counter clockwise, 4 complete weight sets with hanger, 4 pulleys with clamps, two set strings, that allows students to investigate forces in either 3 or 4 forces, durable aluminum table on a large tripod aluminum base with levelling screws.
- Disc diameter: 400 mm
- Set of slotted weights: (50gx9nos + 50g hanger =500g).
- Supplied complete with 4 sets of slotted weight, made of brass.

---

Didac Scientific Pvt. Ltd,  
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom  
Direct Contact Details ✉ [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)  
💻 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)