

Product Code . DS-P-11217

Determination Of Viscosity Of Glycerin



Description

Catalogue No. : "IDAC8303/020"

DESCRIPTION :-

- Viscosity of Glycerine

Contents:

- Digital timer Glass tube
- Ferrite magnet - Nos. 3 'A' base Electromagnet Glycerin - Nos. 5
- Flexible plug lead, 2m black Flexible plug lead, 2m red Steel ball 5.5mm dia.

- Steel ball 13mm dia. Steel ball 19mm dia.

- USB Cable

- Solenoid driver with adaptor 9V

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/