



Product Code . DS-P-11218

Determination Of Density Of Solid

Description



Catalogue No. : "IDAC7304/021"

DESCRIPTION :-

- Density Cube Set with hook (Set of 6) (Brass, Copper, Iron, Aluminum, Wood and Lead)
- Mild Steel Base
- Size : 200X125mm
- With Rubber Feet MS rod (600X 12mm)(Length x Diameter)
- Boss Head Square Bilateral Steel Rod with Hole (Length=200mm) Cotton Thread reel
- Spring balance tubular
- Capacity : 100gm, 250gm
- Beaker
- Material : Plastic
- Capacity : 250ml
- Chemical resistant
- Measuring Cylinder Plastic 250ml
- These autoclavable cylinders have raised prominent graduations for easy reading and have excellent contact clarity.

- The broad base provides good stability.

- Capacity: 250ml.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/