



Product Code . DS-P-11219

Elastic Constant

Description

Catalogue No. : "IDAC1306/022"

DESCRIPTION :-

- For study of Poisson 's ratio etc.,

Salient Features:

- Stainless steal rectangular rods.
- Unyielding wall bracket.

Scope of Supply:

- Stainless Steel rectangular bars each of length 30cm - Nos. 2
- Vernier calliper Screw gauge Stopwatch Thread reel Weighing balance Meter scale
- Experimental wire each L=30cm,dia 1.2mm,1.6mm and 2.5mm - Nos. 3 Wall bracket

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/