



Product Code . DS-P-11221

Newton's Cradle

Description

Catalogue No. : "IDAC7310/024"

DESCRIPTION :-

- It demonstrates that for every action there is an equal and opposite reaction.
- Newton's Third Law of Motion Lift and drop one ball to see one ball pop out from the other side.
- Lift two balls and two pop out.
- 5 steel balls dia. 19mm.
- Mounted on Wooden Base.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/