

Product Code . DS-P-11223

Timing Pro Timing Ball



Description

Catalogue No. : "IDAC9315/026"

DESCRIPTION :-

- The Timing ball is a unique device that records time from moment of release of ball to the moment of impact (on ground or object).
- It can be used to conduct experiments involving acceleration due to gravity (g). Seven Segment FND Display, Rechargeable Battery

Scope Of Supply:

- Timing Ball
- Launching Pad
- Meter Tape 3m

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
🌐 www.didacscientific.co.uk/