

Product Code . DS-P-11225

Boyle's Law Apparatus





Description

Catalogue No. : "IDAC6322/028"

DESCRIPTION :-

- A horizontally mounted acrylic tube having a graduated scale (0-300cm³) along one side.

-
- A gas tap closes the tube.
 - Pressure gauge is present next to gas tap.
 - By changing the position of rubber stopper, the effective length of the tube can be selected.
 - The scale measures volume and pressure by means of pressure gauge.
 - Boyle's law is investigated by changing position of the stopper.
-

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/