

Product Code . DS-P-11226

## Inclined Plane Apparatus



### Description

---

Catalogue No. : "IDAC7323/029"

#### DESCRIPTION :-

- Inclination of plane board can be done and fixed at desired angle.
- The angle may be measured on the angle measurer.
- A pulley is attached at the end, which is used in attaching masses.
- Supplied with roller, pan and two wooden blocks with hook.
- WEIGHT BOXES (PHYSICAL)
  - These are general purpose weights made of nickel plated brass. Provided fractional weights are made of Aluminum.
  - The fraction weights are housed in a separate plastic box.
  - Forceps is provided.
  - Housed in wooden box.
  - Set of weights 1mg to 100gms (no.of weights 9 + fractional weights)
  - Spring balance tubular plastic.
  - Made of transparent acrylic pipe. 500g / 10gm/5n

---

Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom

Direct Contact Details  [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)

 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)