

Product Code . DS-P-11227

Hydraulic Press



Description

Catalogue No. : "IDAC8802/030"

DESCRIPTION :-

- Used to Study The Principle of Hydraulic Press
- It has two graduated syringes situated on a metal frame, having capacities of 50ML and 10ML respectively with pistons.
- Both the syringes are connected through a 3-position tap.
- These syringes may be exposed to atmosphere through this tap, isolated or interconnected as desired.
- Both the pistons are provided with platforms for loading.
- This item is used for demonstrating how fluid pressure is transmitted and mechanical advantage can be achieved.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/