



Product Code . DS-P-11228

Liquid Level Apparatus

Description

Catalogue No. : "IDAC5804/032"

DESCRIPTION :-

- It is used to demonstrate that the level of liquid is same irrespective of size or shape of vessel.
- A common horizontal tube is connected with 4 nos. of vertical tubes (different in shape & size).
- The unit is mounted on a wooden base 230 x 100 x 22 mm, total height 160mm

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/