

**Product Code . DS-P-11229**

## **Spouting Cylinder (Three Hole)**



---

## Description



---

Catalogue No. : "IDAC3809/033"

DESCRIPTION :-

- It is used to demonstrate that pressure increases with the depth.
- It consists of a cylinder with three 4mm dia tube in the same side from top to bottom.
- Dimensions: 40 x 5 cm (Length x Diameter)

---

Didac Scientific Pvt. Ltd,  
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom  
Direct Contact Details  [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)  
 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)