



Product Code . DS-P-11233

Calorimeter

Description

Catalogue No. : "IDAC4827/037"

DESCRIPTION :-

- To determine mechanical equivalent of heat with Joule's calorimeter.

Contents :-

- Joule's calorimeter
- Digital voltmeter: 19.99 VDC
- Digital ammeter: 19.99 ADC
- Battery eliminator: 1. 5 to 12V AC/DC, 3A
- Rheostat: 10 Ω / 3.3A
- Stop watch Thermometer
- Connecting leads red & black 50cm (pair)

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/