



Product Code . DS-P-11234

## Spring Balance Tubular

### Description

---

Catalogue No. : "IDAC5825/038"

DESCRIPTION: -

- Used To Verify Newton's 3rd Law
- Spring Balance Tubular Plastic Set of 6.
- Made of Transparent acrylic pipe.
- Spring mechanism is clearly visible.
- Provided with zero adjustment and two hooks.
- Scaled and calibrated in both grams and newtons.
- Capacity 250gm/2.5N Blue, 500gm/5N
- Green, 1000gm/10N Beige, 2kg/20N Red, 3kg/30N White & 5kg/50N Yellow.

---

Didac Scientific Pvt. Ltd,  
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom  
Direct Contact Details ✉ [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)  
🌐 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)