



Product Code . DS-P-11238

Laser Ray Kit

Description

Catalogue No. : "IDAC5257/042"

DESCRIPTION:

- Laser ray box consists of five independent laser modules with peak wavelength 635nm.
- The laser is very effective in the demonstration of light trace.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/