



**Product Code . DS-P-11241**

## **Optical Bench Hexagonal**

### **Description**

---

Catalogue No. : "IDAC3735/045"

#### **DESCRIPTION :-**

- Length : 1meter
- A high-quality precision bench suitable for advanced work.
- Made of Aluminum alloy with hexagonal section carriage for accessories which can be easily located.
- Accepts the full range of optical accessories with 10 mm diameter mounted rods.
- It is supported at the ends, which ensures accuracy and stability.
- An additional support at the center minimizes the possibility of sagging.
- Supplied complete with 4 fixed saddles.

---

Didac Scientific Pvt. Ltd,  
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom  
Direct Contact Details ✉ [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)  
💻 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)