



Product Code . DS-P-11243

Student's Spectrometer

Description

Catalogue No. : "IDAC5444/46.1"

DESCRIPTION :-

- INTERMEDIATE SPECTROMETER
- It is a robust instrument with high performance.
- A double ended vernier is provided.
- Both telescope and table have a fine adjustment.
- Spectrometer table has lines marked to assist positioning of prism.
- Supplied complete with prism, magnifier, allen key, grating holder & prism holder.
- Scale: Brass scale, Dia 150mm.
- Objective: Achromatic, focal length 178mm, aperture 32mm
- Slit: German silver.
- Reticle: 90° cross on glass.
- Least count: 1 minute.
- Base: Aluminum Casting.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/