



**Product Code . DS-P-11246**

## **Newton's Color Disc**

### **Description**

---

Catalogue No. : "IDAC4575/047"

#### **DESCRIPTION :-**

- It consists of a disc with spectrum colors.
- Diameter of disc 20cm.
- Mounted on stand with turning wheel and handle.

---

Didac Scientific Pvt. Ltd,  
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom  
Direct Contact Details ✉ [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)  
💻 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)