



Product Code . DS-P-11247

Sonometer

Description

Catalogue No. : "IDAC7877/048"

SALIENT FEATURES:

- Pulley with small friction.
- Sharp bridges.
- Tuning fork

Scope of Supply:

- Wooden Sonometer 110 x 10 x 8.5cm³
- Hanger 500g + 6 X 500g slotted weights Rubber Pad
- Screw Gauge
- Tuning forks (512,480,440Hz) with hammer

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/