




Product Code . DS-P-11251

Plane Mirror

Description

Catalogue No. : "IDAC2177/50.3"

DESCRIPTION :-

- Mirror for use in optical experiments.
- SIZE 100X50mm
- Material : Borosilicate Glass
- Shape : Rectangular
- Color : Clear

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/