



**Product Code . DS-C-11259**

## **Boss Head S type**

### **Description**

---

Catalogue No. : "IDAC1255/2.2"

DESCRIPTION :-

- Boss Head S Type
- Double ended, made of Aluminum for holding square/round rod perpendicular to each other.
- Clamp opening 16mm.
- High Efficiency
- Anti-corrosive
- 'T' thumb screws made of chromed Iron.

---

Didac Scientific Pvt. Ltd,  
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom  
Direct Contact Details ✉ [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)  
💻 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)