

Product Code . DS-C-11261

Bunsen Burner





Description

Catalogue No. : "IDAC7452/04"

DESCRIPTION :-

- A Bunsen Burner is a mechanical apparatus that is connected to a flammable gas source.
- There is a knob to adjust the amount of gas flow and a rotating collar that controls airflow.
- These both would be adjusted to get an ideal flame for heating purposes.
- Nickel plated burner tube of 12mm dia with air regulator and a riffled connector.
- Mounted on a painted cast iron base.
- Only the air supply port can be adjusted.

-
- without stopcock

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/