



Product Code . DS-C-11265

Porous Plate

Description

Catalogue No. : "IDAC2156/8"

DESCRIPTION :-

- Porcelain cavity tiles
- Chemical resistant
- Made of Porcelain.
- Glazed, cavity used for spot test.
- It can hold liquid or precipitate in each well at the time of reactions.
- No of Depression 6

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/