



Didac Scientific



Product Code . DS-C-11270

Small Thin Walled Test Tube



Description

Catalogue No. : "IDAC9835/13"

DESCRIPTION :-

- Test tube without rim
- Size:10x75mm
- Size : Small
- Type : Thin walled
- Material: borosilicate glass
- Color: clear
- Chemical resistance
- Heat resistance?

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/