



Product Code . DS-C-11278

Measuring Cylinder



Description

Catalogue No. : "IDAC2144/18"

DESCRIPTION :-

- Material: Polypropylene

-
- Capacity: 50ml
 - The Measuring Cylinders, moulded in Polypropylene & TPX, are single piece cylinders.
 - These autoclavable cylinders have raised prominent graduations for easy reading and have excellent contact clarity.
 - The broad base provides good stability.
-

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/