



**Product Code . DS-C-11279**

## **Test Tubes Stand**

### **Description**

---

Catalogue No. : "IDAC8455/20"

#### **DESCRIPTION :-**

- Test tube stand plastic 6 holes 6 pegs
- This autoclavable Stand is perhaps the most commonly used Test Tube Stand in various laboratories.
- For better visibility, this 6-place rack is end supported by 2 columns.
- The top plate has 2 holes for holding 25 mm dia test tubes, four holes to hold 16 mm dia tubes whereas the base has corresponding hemispherical wells.
- It also has six vertical pins for drying tubes. Another option available in this section is a test tube stand in which all the holes in the top
- plates can hold 25 mm dia tubes.

---

Didac Scientific Pvt. Ltd,  
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom  
Direct Contact Details ✉ [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)  
💻 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)