



Product Code . DS-C-11281

Weighing Bottle

Description

Catalogue No. : "IDAC4112/22"

DESCRIPTION :-

- Size: 60x25mm.
- Material: borosilicate glass.
- With interchangeable glass stopper.
- Chemical resistance
- Heat resistance
- Color: transparent

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/