

Product Code . DS-C-11286

Standard Flask





Description

Catalogue No. : "IDAC4122/27"

DESCRIPTION :-

- Capacity: 250ml
- Material: polypropylene

-
- Strong, unbreakable and have good contact clarity.
 - These autoclavable flasks are provided with screw caps to prevent leakage.
-

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/