

Product Code . DS-C-11288

## Burette with Stand



### Description

---

Catalogue No. : "IDAC8541/29"

Features:

- Best quality raw material
- Long life
- Most appreciated

---

## DESCRIPTION :-

- Retort stand 225 x 150mm side hole.
- These non-corrosive, moulded Retort Stands provide excellent replacement for the metal ones.
- Even the metallic rod has a polypropylene covering making it non-corrosive.
- The rod measuring 12.5 mm x 750 mm can be screwed into the heavy base giving way to a perfectly stable Retort Stand.

---

Didac Scientific Pvt. Ltd,  
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom  
Direct Contact Details ✉ [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)  
💻 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)