



Product Code . DS-C-11291

Measuring Jars

Description

Catalogue No. : "IDAC3256/30"

DESCRIPTION :-

- Capacity: 500ml
- Material: Polypropylene
- These handy Jugs, moulded in Polypropylene, are clear, autoclavable and have good chemical resistance.
- These jugs are provided with a handle for easy & convenient lifting and have raised graduations for easy reading.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/