

Product Code . DS-C-11294

## Funnel



### Description



---

Catalogue No. : "IDAC6589/33"

DESCRIPTION :-

- Analytical funnel polypropylene
- Designed specifically for analytical chemistry, the body of the funnels is 600.
- Moulded in PP, these funnels are provided with internal & external ribs to prevent air blockage.
- Diameter: 75mm.
- Inside Dia: 75mm
- Stem Length: 77mm

---

Direct Contact Details  [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)  
 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)