



Product Code . DS-C-11295

Glass Rod

Description

Catalogue No. : "IDAC7422/34"

DESCRIPTION :-

- Glass rod
- Material: borosilicate glass
- Length dia: 6 X 300mm
- Color: clear
- Chemical resistance
- Heat resistance

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/