

Product Code . DS-C-11298

## Iron Stand with Clamp



### Description

---

Catalogue No. : "IDAC1027/37"

#### DESCRIPTION:-



- A Retort Stand is a piece of scientific equipment to which clamps can be attached to hold other pieces of equipment and glassware for instance, burettes, test tubes and flasks.
- Made of cast iron powder coated mild steel rod chrome plated.
- Retort Stand Base (without Rod) 200x125mm.

- 
- Material: mild steel sheet Powder coated finish.
  - Supplied without rod.

#### FEATURES: -

- Best quality raw material
- Long life
- Most appreciated

---

Didac Scientific Pvt. Ltd,  
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom  
Direct Contact Details  [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)  
 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)