

Product Code . DS-C-11302

Crystallizing Dish



Description

Catalogue No. : "IDAC5896/39"

DESCRIPTION :-

- Color: clear
- Chemical resistance
- Heat resistance
- Material: borosilicate glass
- Size: dia x height 100x50mm

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/