

Product Code . DS-C-11309

Rubber Stopper



Description

Catalogue No. : "IDAC4596/41.5"

DESCRIPTION :-

- Without hole no.7
- Top diameter 26.0mm
- Bottom diameter 20.0mm
- Height 27.5mm.
- Material : Rubber
- Chemical resistant

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/