



**Product Code . DS-C-11327**

## **Glass Beaker**

### **Description**

---

Catalogue No. : "IDAC5827/011"

#### **DESCRIPTION :-**

- Material: polypropylene
- Capacity: 250ml
- Beakers, molding Polypropylene, have excellent clarity and very good chemical resistance.
- These autoclavable Beakers have prominently raised graduations on the outer surface.
- Pouring a solution becomes very convenient because of its tapered spout.

---

Didac Scientific Pvt. Ltd,  
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom  
Direct Contact Details ✉ [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)  
🌐 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)