

Product Code . DS-C-11329

Blue Litmus Papers



Description

Catalogue No. : "IDAC7431/13"

DESCRIPTION :-

- The blue litmus will turn red or pink when it comes into contact with an acidic solution below pH 4.5.
- Between 4.5 and 8.3 the paper will be purple.
- Above 8.3 it will stay blue.
- Supplied in plastic vial, about 100 paper strips/vial
- Blue – color changes from Blue to red when dissolved in base solution.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk

