



**Product Code . DS-C-11331**

## **Bunsen Burner**

### **Description**

---

Catalogue No. : "IDAC9840/15"

#### **DESCRIPTION :-**

- Burner bunsen without stopcock
- Common piece of laboratory equipment that produces a single open gas flame, which is used for heating, sterilization, and combustion.
- Nickel plated burner tube of 12mm dia with air regulator and a riffled connector.
- Mounted on a painted cast iron base.
- Only the air supply port can be adjusted. without stopcock

---

Didac Scientific Pvt. Ltd,  
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom  
Direct Contact Details ✉ [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)  
💻 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)

