




Product Code . DS-C-11332

Wire Gauze

Description

Catalogue No. : "IDAC1369/16"

DESCRIPTION :-

- Size: 15 x 15cm
- Shape: Square
- Ceramic Circle: Provided in center
- Wire Edges: Rolled Smooth
- Filled centre for use on tripods.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/