



Product Code . DS-C-11333

Tripod Stand



Description

Catalogue No. : "IDAC2847/17"

DESCRIPTION :-

- Tripod stand triangular (cast top).
- Height x side of triangular 15x12.5cm
- Material: Cast iron top and mild steel legs painted.
- Zinc coated legs are bent outwards for stability.
- A tripod is a portable three-legged frame, used as a platform for supporting the weight and maintaining the stability of some other object.

-
- A tripod provides stability against downward forces and horizontal forces and movements about horizontal axes.
 - A laboratory tripod is a piece of three-legged equipment commonly used to conduct experiments in laboratories.
 - It is used as a platform to hold and support glassware, such as beakers and flasks, during experiments
-

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/