

Product Code . DS-C-11337

## Knife



### Description

---

Catalogue No. : "IDAC6782/21"

DESCRIPTION :-

- It is made of Stainless-Steel
- Having size no. 4
- Chemical resistance
- It can be used with blade with blade no. 20, 21, 22, 23, 24, 25. (Only handle without blade)

---

Didac Scientific Pvt. Ltd,  
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom  
Direct Contact Details ✉ [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)  
🌐 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)